TISSUE PACKAGE INSERT QC-605-F-107 V 2.0

DESCRIPTION
 DONATED HUMAN TISSUE. Tissue grafts are recovered from deceased human donors. The procedure executed to manufacture the graft, including recovery, donor screening, testing, processing, packaging, labeling, storage, and distribution were performed in compliance with all applicable local, state, and federal regulations, including the use of Tissue Banks Standards for Tissue Banking and the applicable European Union Directives for human tissue. Tissue is processed in a clean room environment, following rigorous quality assurance standards. Tissue labeled as **STERILE R** has been terminally sterilized to a SAL of 10⁻⁶ (Sterility Assurance Level) using Gamma irradiation (Cobalt 60).

SCREENING AND TESTING
 The donor has been determined to be eligible by a Community Tissue Services Medical Director at 349 S. Main St., Dayton, OH 45402 based on the results of screening and testing. Screening includes a review of medical and social history, hospital records, infectious disease screening, autopsy report (if performed), and all other relevant information. All donors are screened for HIV-1, HIV-2, HBsAg, Anti-HBc, Anti-HCV, HCV NAT, HBV NAT, HCV NAT, and Syphilis (RPR, VDRL, RIA, and TPP). All donors are screened for NAT and syphilis. FDA-approved test kits are used when available. Additional tests, including but not limited to HTLV I/II, may have been performed and were found to be acceptable for transplantation. Communicable disease testing has been performed by a laboratory registered with the FDA to perform donor testing and certified to perform such testing on human specimens in accordance with the Clinical Laboratory Improvement Amendments of 1988 (CLIA) and 42 CFR part 493, or that has met equivalent requirements as determined by the Centers for Medicare and Medicaid Services (CMS).

STORAGE
 Freeze-Dried tissue must be stored at ambient temperature. It is the responsibility of the tissue dispensing service, tissue distribution intermediary, and/or end user clinician to maintain tissue intended for transplantation in appropriate storage conditions prior to further distribution or transplant. Tissue may not be stored at liquid nitrogen (LN) vapor phase or LN liquid temperatures.

AREAS OF APPLICATION
 Human bone screw grafts are used for implantation in children and adults.
 Shark Screw™ ACL: Used for fixation and reconstruction of tendons and ligaments (e.g. ACL, PCL).

PRODUCT SIZES AND APPLICATION INSTRUMENT COLOR CODE
 Ø 3.5 mm cut; hexagonal head
 Ø 3.5 mm, 35 mm length; blue color code
 Ø 4.0 mm, 35 mm length; yellow color code
 Ø 4.5 mm, 35 mm length; red color code
 Ø 5.0 mm, 35 mm length; black color code
 Shark Screw™ ACL: hexagonal head
 Ø 7.0 mm, 21 mm length; cannulated; light green color code
 Ø 8.0 mm, 21 mm length; cannulated; light green color code
 Ø 9.0 mm, 21 mm length; cannulated; light green color code

The appropriate Shark Screw™ variant and size must be selected by the surgeon and depends on the indication, the anatomical region, and the associated mechanical load as well as the size of the host bone.
 The Shark Screw™ grafts can only be applied with the use of the corresponding surgical instrumentation. These are available as separate products. Instrument accessibility and readiness for use must be ensured.

CONTRAINDICATIONS

- Insufficient bone substance to anchor the graft.
 - Necrotic host site.
- PRECONDITIONS FOR RESTRICTED INDICATION**
- Insufficient bone quality, bone formation disorders (e.g., Osteoporosis). The bone quality must be checked in advance.
 - Application in the host site with poor blood supply or in the infected host site (due to a poorer healing rate).
 - Circulatory disorders slowing down the healing process.
 - Acute or chronic infections.
 - Circumstances which prevent the patient from restricting their activities accordingly or from following the doctor's instructions during the healing phase.
 - The graft must not impair the growth plate.

PREPARATION

- The grafts are in a double layer peel package: peel outer package down and aseptically deliver inner package to the sterile field or sterile team member.
- Open the second peel pack right before application and apply within 30 min after opening.
- The graft is used **without prior rehydration**. The graft swells in the host bed by approximately 2%.

APPLICATION OF SHARK SCREW™ cut

- The graft keeps the compression stable after being secured in.
- The corresponding, color-coded instruments must be used for the preparation of the insertion channel. After drilling the hole the thread is cut with the matching thread cutter. Deviation from the original drill hole axis must be avoided. The finished screw channel must be rinsed to remove any remaining bone debris. The screw is inserted with the appropriate insertion instrument. The optimal graft position and depth is controlled by imaging if necessary.

QC-605-F-107 V 2.0

- After insertion, the screw graft cut is chlorinated to the bone head.
- The hexagonal head of the Shark Screw™ cut needs to be cut off after insertion. The hexagonal head must not be driven into the bone (danger of bursting).

APPLICATION OF SHARK SCREW™ ACL

- The grafts are designed as interference screws with a hexagonal screw head that can be countersunk underneath the bone surface.
- Preparation of the bone surface: The correct screw graft, and the thickness of the tendon to be fixed, the appropriate diameter of the drill hole is to be determined by the surgeon. The drill hole must be at least the diameter of the used screw graft.
- No thread is cut.
- 1.2 mm wire may be used for guiding while screwing in the screw graft.
- The tendon/ligament to be fixed is pulled through the drill hole.
- The graft is inserted beneath the bone surface with the insertion instrument marked in light green. The tendon/ligament is fixed by the resulting contact pressure.

GENERAL NOTICE ON SHARK SCREW™ APPLICATION:

- The risk that the bone screw graft breaks during insertion is very low. If the graft should break, the inserted portion may remain in the bone. A new hole may be drilled next to the graft or through the graft (cut) and a new graft can be inserted.
- It is necessary to immobilize the affected area after surgery to ensure healing of the graft. The duration of immobilization depends on various parameters such as localization of the graft, type of operation, age of the patient, bone quality and bone wound area. The exact duration of immobilization is at the discretion of the treating physician. The immobilization phase is followed by the mobilization phase for the patient.
- Before using the Shark Screw™ graft, it is recommended that the surgeon be acquainted with the product's specific surgical methods, information and operation videos are available on www.surgebright.com.

WARNINGS AND PRECAUTIONS

- Intended for use in one patient, on a single occasion only.
- Do not use if package integrity has been compromised. Once the user breaks the container seal, the tissue grafts must be transplanted or discarded.
- Tissue may not be sterilized or re-sterilized.
- This tissue is intended for use by qualified healthcare specialists such as physicians, dentists, or podiatrists.
- Although this tissue has been tested and screened for human pathogens, and processed under aseptic conditions, human derived tissue may still transmit infectious agents.
- Community Tissue Services (US only) or SurgeBright (EU only)
- Tissue has been processed with betadine and/or Polymyxin B and traces may remain.

TISSUE TRACKING

Recipient records must be maintained for the purpose of tracing tissue post-transplantation. Tissue traceability records including the SEC (applicable to EU only) must be maintained by the organization responsible for human application (ORHA). Information on Community Tissue Services (applicable to grafts implanted in the US only) Federal Regulations (21 CFR 1271) and Joint Commission Standards require exact tracking of the tissue. It is the responsibility of the end user to provide this information, which enables Community Tissue Services to maintain records for the purpose of tracing the tissue post-transplant.
 Community Tissue Services is accredited by the American Association of Tissue Banks. Community Tissue Services –Center for Tissue, Innovation and Research is ISO 13485 certified. Health Canada Registration: 100076.
 SurgeBright and Community Tissue Services make no claims concerning the biological or biomechanical properties of the tissue. Community Tissue Services declares all liability and responsibility for any misuse of tissue intended for clinical application.
 Please contact SurgeBright GmbH or Community Tissue Services should you require further information.

Processed by:
 Community Tissue Services – CTR
 Center for Tissue, Innovation and Research
 Manufacturing and Distribution Center
 2900 College Drive
 Kettering, Ohio 45420
 (800) 684-7783
 Fax (937) 461-4237
 Fax (937) 222-0228

SYMBOL GLOSSARY	
	Do not reuse
	Do not resterilize
	Sterilized using irradiation
	Consult instructions for use

Released and Distributed by:
 CTR
 Center for Tissue, Innovation and Research
 Manufacturing and Distribution Center
 2900 College Drive
 Kettering, Ohio 45420
 (800) 684-7783
 Fax (937) 222-0228

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Hergestellt durch:
 Community Tissue Services – CTR
 Manufacturing and Distribution Center
 2900 College Drive
 Kettering, Ohio 45420, USA
 (800) 684-7783
 Fax (937) 222-0228

Freigegeben durch:
 Community Tissue Services – CTR
 Manufacturing and Distribution Center
 2900 College Drive
 Kettering, Ohio 45420, USA
 (800) 684-7783
 Fax (937) 461-4237

SYMBOL GLOSSARY	
	Do not reuse
	Do not resterilize
	Sterilized using irradiation
	Consult instructions for use

How to return this form	
Email	tissueusage@patienttracking.care
Fax	937-222-2538
Mail	Community Tissue Services Attn: Tissue Usage 2900 College Dr. Kettering, OH 45420

Allograft Tracking Form

FDA Regulations and Joint Commission Standards require tissue tracking systems in all hospitals using allograft tissue for transplantation. In order to comply with these requirements, please complete ALL fields on this form.

Date of Surgery: _____

Patient's Medical Record Number or Date of Birth: _____

Community Tissue Services does not consider the information requested on this form to be protected health information (PHI), as defined under the HIPAA regulations. Information considered to be PHI by the originator should not be released to Community Tissue Services.

**Place peel-off label for up to 4 allografts or write tissue ID# in the spaces provided.
One patient, one procedure per tracking form.**

Allograft Tissue ID# _____ Place Peel-Off Label Here
--

Allograft Tissue ID# _____ Place Peel-Off Label Here
--

Allograft Tissue ID# _____ Place Peel-Off Label Here
--

Allograft Tissue ID# _____ Place Peel-Off Label Here
--

**If any questions, problems, or adverse reactions occur,
contact 1-800-684-7783 or 1-937-222-0228.**